



EXPLORE
Minimally Fractional RF Technology



#1206 SAMSUNG IT Valley, 197-5, Guro-dong, Guro-gu, Seoul, Korea. Tel: 82-2-6925-4459(Rep.) Fax: 82-2-6925-4460
E-mail: danalaser92@gmail.com Website: www.i-dana.com



www.i-dana.com

EXPLORE

Minimally Fractional RF Technology



Fine Wrinkle · Pore · Acne Scar · Skin Elasticity · Neck Wrinkle



Specification

Type	Electro-surgical Device
Output Power Level	1Level ~ 120Level
Needle Length	0.1mm ~ 4 mm
Repeat Time	0.1 sec ~ 3.0 sec
Needle type	Invasive / Non-invasive
Electrical	100~240VAC, 50Hz/60Hz
Dimension	330(W)×490(D)×1175(H)mm
Weight	23kg



Minimally Fractional RF Technology

High Treatment Effect

Minimize pain

Less Down Time

No side effect

Cover all ages




The harmony of 2 way Hand-piece



Effect of EXPLORE

Collagen remodeling in Dermal through Wound Healing Process



Minimally Fractional RF Technology



- Model : EXPLORE
- Manufacturer : DANA
- Supplier : DANA
- Application : Face & Neck Lifting, Wrinkle, Hyperhidrosis, Whitening, Scar, Pore, Acne